Simplewell 昇微

CO2 incubator

Simplewell Technology Co., Ltd

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Email: sales01@simplewell.com.cn

Contents

- 01.** Production description
- 02.** Product Innovation Features
- 03.** Advanced technical indicators
- 04.** Customer promotion

01
Part

Production description

1.1 Product application

Scope of application

Suitable for cell, tissue, bacterial culture, immunology, oncology, genetics and bioengineering, and widely used in research and production in microbiology, medicine, pharmaceuticals, environmental protection, food, animal husbandry etc. Its principle is to use artificial methods to create an artificial environment for the growth and reproduction of microorganisms, cells, bacteria, and viruses in an incubator, such as controlling a certain temperature, humidity, gas, etc., and is mainly used for tissue culture and the cultivation of some special microorganisms.

Appearance


Panoramic glass interior door

Coated inner cabin


1.2 Operation interface

Main interface

Setting interface

Run interface

1.3 Heating system


Heating tube

Customized heating tube, Long life and reliable quality

PT100 temperature sensor

Overtemperature protection

Temperature level one protection

Solid state relay

Brand: Carlo gavazzi Swiss
To control heating tube interruption

Gas protection switch

Brand: Rainbow Kirea
Temperature level two protection

- 1.Fin-tube heater
- 2.With dry heating prevention device
- 3.Heating control: adopt PWM pulse intelligent width adjustment control technology

1.4 CO₂ system

Co2 sensor

high temperature resistance and has no impact on performance when cleaning at 180°C.

Float flowmeter

High adjustment accuracy and long service life

Gas pipe

Soleniod valve

Brand Airtac, reliable quality and durable

CO₂ system: Adopt high-temperature resistant CO₂ sensor, no need to disassemble for high-temperature cleaning

1.5 Electric control

Display

Self-developed controller

USB

Socket

Standard plugs in line with national standards

Flame retardant wire

Motor

Customized motor, safe and reliable

Leakage switch

Schneider

Power wire

Standard plugs in line with national standards

The control parts such as air switch, contactor, relay, touch screen, frequency converter, etc. adopt international famous brands

1.6 Structure component

Balance goblet

。

Door hinge

Humidification
water tray

Double sided
protective coil


Adjustable height sample
rack

Cooling tube

Can move up and down to
facilitate the movement of the
water tray

Throttle

Standard throttle, easy to replace and
highly replaceable

Tempered glass

1.7 Standards compliant

- GB/T 28851-2012 Specifications for refrigerated incubators
- YY 1621-2018 Medical carbon dioxide incubator
- DIN 12880

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Appearance, structure, technical performance indicators of the biochemical incubator specified in GB/T 28851-2012

Appearance, structure, temperature and humidity control, carbon dioxide concentration etc specified in YY 1621-2018

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1.8 Process advantage

1. Pipeline welding process: high-quality copper tube nitrogen shielded welding method is adopted, which avoids the damage to the compressor caused by the oxide impurities on the inner wall of the copper tube entering the refrigeration system caused by the traditional welding method.

3. Pipeline protection measures: The pipeline of the refrigeration system adopts the method of adding anti-vibration hose and C-shaped elbow to avoid copper pipes and cracks caused by vibration and temperature changes.

5. When the equipment is running, detect the circuit temperature of the power distribution cabinet.

2. Damping measures: Install damping springs and anti-vibration soft rubber pads at the bottom of the compressor and pipeline to reduce vibration.


4. Noise control: The condenser adopts the German Marl low-speed high-air volume condensing fan, and installs wave-shaped sound-absorbing sponge around the refrigeration unit to achieve lower noise.

1.9 Manufacturing process

Vibration test of key components
simulated road conditions

Parts such as evaporator, motor and water tank are
subjected to vibration tests before installation

Vibration test of the whole machine
simulating bumpy road conditions

Vibration test for small equipment
before shipment

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
[illegible]

Confirmation of the production process: After the production of the equipment begins, the person in charge of each link will carefully fill in the confirmation form to correct the problems in the production process in a timely manner. At the same time, trace the source, optimize the production process, improve production efficiency, and ensure the quality of each equipment produced

1.11 Manufacturing process and requirements

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Flame-retardant storage board

The picture shows the performance test report of the storage board in terms of flame retardancy, compressive strength, and bending bearing capacity

1.12 Manufacturing process


Adopt flame retardant wires, the picture shows the wire flame retardant certification

02
Part

Features

2.1 Product innovation features

Simplewell CB series internal air duct circulation


Simplewell adopts a pollution-free jacketed air supply duct structure, which effectively improves the temperature uniformity of the equipment.

Jacketed air duct and circulation direction

Experimental cabin

Boss formed by stamping and drawing

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Advanced technical indicator

3.1 Advanced related technologies

- The external flanging of the Simplewell instrument is stretched by one-piece mold, with natural transition
- The inner wall of the SimpleWell equipment is stamped and stretched to form a boss, which is beautiful and easy to use. It not only improves the strength of the equipment, but also facilitates the placement of opposite shelves, reduces welding points, screw fixing points and other unnecessary connectors, and enhances the Background of sterile cleanliness of equipment.
- Simplewell adopts a large field of view tempered glass window, which effectively enlarges the field of view to observe the situation, and can quickly find the situation of the medicine inside.
- Simplewell uses special materials to coat the inside of the equipment and shelves to facilitate high-temperature cleaning and internal cleaning. There will be no residual water vapor hanging on the wall, effectively avoiding internal environmental pollution and improving cleanliness.
- 180°C high temperature sterilization

Integrated flanging

Integrated stamping boss

Large view glass window

Special internal coating


04

Part

Customer promotion

4 Customer case

Simplewell can customize production to meet the special needs of different customers

Simplewell 昇微

Thanks for watching

Simplewell Technology Co., Ltd

Address: No.221, Shuixin Road, Dalang Town, Dongguan City

Tel: 0769-88887909 Fax: 0769-88885229

Website: www.simplewell.com.cn

Email: sales01@simplewell.com

Team

Persistence

Cooperation

Honor